

Retrotec DU200 DucTester

Fan Press. Pa	Open High CFM	Mid CFM	Low CFM
10	100		
15	123		
20	142	30	6
25	160	34	7
30	175	37	8
35	190	41	9
40	203	44	10
45	216	46	10
50	228	49	11
55	239	52	12
60	250	54	12
65	260	57	13
70	270	59	14
75	280	61	14
80	289	63	15
85	299	66	15
90	307	68	16
95	316	70	16
100	324	72	17
105	333	74	17
110	341	75	18
115	349	77	18
120	356	79	19
125	364	81	19
130	371	83	20
135	378	84	20
140	385	86	21
145	392	88	21
150	399	89	22
155	406	91	22
160	413	92	23
165	419	94	23
170	426	95	23
175	432	97	24
180	438	98	24
185	444	100	25
190	451	101	25
195	457	103	25
200	463	104	26
205	468	106	26
210	474	107	27
215	480	108	27
220	486	110	27

Fan Press. Pa	Open High CFM	Mid CFM	Low CFM
230	497	112	28
235	502	114	28
240	508	115	29
245	513	116	29
250	518	118	30
255	524	119	30
260	529	120	30
265	534	121	31
270	539	123	31
275	544	124	31
280	549	125	32
285	554	126	32
290	559	127	32
295	564	129	33
300	569	130	33
320	588	134	34
340	607	139	36
360	625	143	37
380	642	148	38
400	659	152	39
420	676	156	41
440	692	160	42
460	708	164	43
480	724	167	44
500	739	171	45
520	754	175	46
540	769	178	47
560	783	182	48
580	797	185	50
600	811	189	51
620	825	192	52
640	839	196	53
660	852	199	54
680	865	202	55
700	878	205	56
720	891	209	57
740	903	212	58
760	916	215	58
780	928	218	59
800	940	221	60
820	952	224	61
840	964	227	62
860	975	230	63